


# COMPACT CHEST SYSTEM



CCWS Wall Stand




CCTS Tube Stand


Compact Installation in:  
Doctor's Office  
Small Clinics  
Chiropractors  
Orthopedics



Patriarch



## GTR Compact Chest System



### System Specifications

Wall Stand:	Vertical Wall Stand, Electromagnetic Locks
CCWS	17" Grid Cabinet/Tray, GRID (103LPI, 10:1) Optional 36"
	Optional Reciprocating Bucky and AEC
X-Ray Tube Support:	Floor / Wall Single Rail Tube Stand Electromagnetic Locks on all movements
CCTS	40" to 72" focal film distance
	Floor Rail Length: 53"
	Arm Length: 24"
	Lowest Focal Spot: 12"
	Highest Focal Spot: 73"
	Tube Rotation: +/- 180 degrees
X-Ray Tube:	
	125 kVp 1.0/2.0 Focal Spots, 140 kHu, 3" Anode with Trunion
	125 kVp 1.0/2.0 Focal Spots, 150 kHu, 3" Anode with Trunion
Collimator:	Manual Collimator for Optimum Protection and Easy Use



GTR Patriarch  
High Frequency X-ray Generator

510 Elk St.  
Gassaway, WV 26624

GTR® Labs, Inc.

[gtr@gtrllc.com](mailto:gtr@gtrllc.com)  
[www.gtrllc.com](http://www.gtrllc.com)

Phone: 1-888-871-9232  
304-364-2211  
Fax: 304-364-2212