


GTR® VET RAD SYSTEM




GTR® offers you the latest in complete radiographic x-ray systems for veterinary applications. Utilizing the GTR® Intelligen® High Frequency X-Ray Generator with the VET APR Console, the VET RAD SYSTEM allows the maximum flexibility both for your budget and your application.



Patriarch Generator



VTT Integrated Table/Tube Stand



VET APR Console



From the basic system pictured on the front of this data sheet to the feature-enhanced system pictured here, GTR® can offer the right combination of equipment to meet your needs.

We have the basic table/tube stand combination, popular with small animal veterinarians, available with a 4-way float top and a variable SID.

For those more demanding applications where 360 degree access is required and floor space must be clear, we have the system which utilizes the overhead tube support and elevating table with 4-way float top.



VOT Overhead Tube Support
Handlebar APR Console



VTE Elevating Table

High Frequency X-Ray Generators:

Fully Tested with a 5 year Warranty

Patriarch 37.5 kW High Frequency X-Ray Generator (400mA/125kVp)
APR Console Mounted in Handlebars
2 Position Footswitch for prep & expose

Tables: *All Tables factory assembled with grid cabinet, stainless steel tray and 103 line, 10:1 grid, 15" Longitudinal Travel of grid cabinet*

VTT Integrated Table/Tube Stand
Table is standard with a fixed top 60"x28" Height 29.5"
Table Base: 19" x 39.5"

Option: 4-Way Float Table Top Table Top Longitudinal Travel: 27" Table Top Transverse Travel: +/- 4" Tube Stand to Table Edge, Fully Extended: 43.5"	Treadle Switch on Table Top Locks Electromagnetic Locks on all movements
---	---

VTE Elevating Table with 4-Way Float Top 60"x28" Height 21.5" to 33.5"

X-Ray Tube Supports:

VTT Integrated Table/Tube Stand
Tube Column Height: 75"
Tube Stand Longitudinal Travel: 40"
Tube Stand Vertical Travel (Variable SID): 29"
Tube Rotation: +/- 90 degrees
Electromagnetic Locks on all movements

VOT Overhead Tube Support with longitudinal travel, SID fixed
APR Console mounted in handlebars

X-Ray Tubes:

1.0/2.0 focal spots, 140 kHu
0.6/1.5 focal spots, 200 kHu

Collimators:

Manual Collimator for optimum protection and easy use

Your Dealer:

GTR® Labs, Inc.
510 Elk St.
Gassaway, WV 26624
Phone: 304-364-2211
Fax: 304-364-2212

gtr@gtrllc.com
www.gtrllc.com